Estimado Cliente, Dear Customer,

Gracias por habernos escogido. Thank you for choosing us.

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3. Material necesario / Required material

1	Pipeta 30 ml 2 aforos Pipette 30 ml 2 marks	ref. 3019132
1	Pipeta 10 ml graduada Pipette 10 ml graduated	ref. 3019017
1	Probeta 25 ml Graduated cylinder 25 ml	ref. 3018002
1	Soporte bureta Burette holder	ref. 1010007
1	Pinza bureta Burette clamp	ref. 63204115

1	Bureta 10ml Burette 10ml	ref. 3019063
1	Vaso 50ml Beaker 50ml	ref. 3026162
1	Contador de tiempo-Timer Timer-Counter	ref. 30880081
1	Indicador flujo Flow indicator	ref. 45000539
1	Soporte matraces PR Flasks holder PR	ref. 1010025

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The official method is to convert sulphites in the wine into undissociated sulphur dioxide by acidifying the wine sample with phosphoric acid to a pH of about 1. The volatile sulphur dioxide is liberated from the sample with the help of an air stream bubbling through the wine sample. When this takes place the 'free' sulphur is liberated. The 'bound' sulphur is not dissociated at room temperature but with gentle boiling dissociates into sulphur dioxide, which is again liberated by the air stream.

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04. Cerrar los grifos de cada matraz reactor. Ensure the faucet of each reactor flask is closed.

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of the wine sample and immediately close with the screw cap

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Switch on the pump and aerate with 1 l/min of air for at least 20 minutes. Check that all flasks receive the same amount of air. Otherwise, try to reposition them or refer to the Observations section.

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Once the required time has elapsed, switch off the pump and unclamp the pear shaped flask. Titrate it with 0.01M Sodium Hydroxide back to the olive green colour. Record the spent volume of the burette and multiply it by 10.7 to obtain the free SO₂ expressed in mg/l.

NOTA / NOTE: El factor será 16 si la concentración de la sosa es exactamente 0.01 M y se usan 20 ml de vino. The factor is only 16 if the concentration of the sodium hydroxide is exactly 0.01M and 20ml of wine are used.

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Proceed identically to free sulphur but with heating of the sample by pressing the switch associated with each distiller. You can also proceed after a free SO₂ analysis. In that case, after the titration of the Free SO₂ replace the flask shape pear in its position and turn on the heating element and the air pump for 20 minutes more. This will be the distillation of the Combined SO₂. Once past time remove the flask and titrate again. Total SO₂ will be the sum of the Free plus the Combined SO₂.

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Ensure the same bubbling (aeration) occurs in each SO₂ bubbler flask (pear shaped flask).

Factorizar periódicamente la sosa y corregir el factor multiplicador. Periodically factorise the sodium hydroxide and correct the multiplier factor.

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When the pump is turned on, make sure that bubbles are produced in all flasks at the same speed. If not, check the upper and lower air outlet connections.

Si, por alguna razón, una de las unidades deja de estar operativa, haga un by-pass con el tubo suministrado de la unidad anterior a la inoperativa a la siguiente. If, for some reason, one of the units is no longer operative, make a by-pass with the tube supplied from the previous unit to the inoperative to the next one.

Tenga cuidado al agregar el ácido orto-fosfórico al reactor. En caso de contacto con la pie neutralícelo con bicarbonato sódico primero y luego lávese las manos con suficiente agua. Be careful when adding orthophosphoric acid to the reactor. In case of split or contact with the skin neutralize it with sodium bicarbonate first and later wash your hands with enough water.

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All instruments have a warranty of 2 years against manufacturing defects except misuse.

La garantía no cubre fungibles tales como resistencias o roturas de vidrio excepto que este tenga defectos de fabricación. Warranty does not cover consumables such as heating elements or breaks of glass except that this have manufacturing defects.

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follow the instructions of the equipment or contact us if any doubts are risen.

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Reclamations will not be accepted after 15 days.

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No returns will be accepted if the instrument is not received in perfect conditions (without stains, scratches, dirtyness,...) and in its original package.

Las devoluciones aceptadas tendrán un demérito del 20%. Accepted returns will have a 20% loss in price.

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Any reparation has a warranty of 6 months.

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Instruments to be repaired will be received at prepaid shipment, well packed and without glass components thus not accepting shipment to be paid. In that case packet will be returned.

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All products are tested and shipped well packed from factory. If the product arrives broken please complain immediately to the courier. GAB Sistemática Analítica S.L. is not responsible for possible damage done by shipment.

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All instruments are tested before being shipped so some might contain traces of water.

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