



# Reducing matters kit (AAZ)

ref. 1002000

The simplest analytic for the determination of sugars in wines, musts, ciders, beers, ... by Rebelein method.

## Characteristics

- Fast
- Easy titration
- Direct result with inverted GAB burette
- Suitable for tirage (2nd fermentation)
- Matrix independence
- Up to 50 determinations
- Requires GAB sugars burette
- Requires heater GAB
- Rebelein method



## Pack composition

- |                              |                |
|------------------------------|----------------|
| · 1 x 500 mL Acid            | ref. 1002001   |
| · 1 x 250 mL Alkaline sol.   | ref. 1002002   |
| · 1 x 500 mL Copper sol.     | ref. 1002003   |
| · 1 x 500 mL Iodine          | ref. 1002004   |
| · 1 x 500 mL Starch          | ref. 1002005   |
| · 1 x 1000 mL Titration sol. | ref. 1002006   |
| · 1 pumice stone tube        | (ref. 1015001) |
| · 1 instructions manual      |                |

## Required material

- |                            |                  |
|----------------------------|------------------|
| · 1 electric heater GAB    | ref. 1010005     |
| · 1 GABA AZ burette        | ref. 3019033     |
| · 1 watch glass 60 mm      | ref. 3026003     |
| · 1 clamp test tubes       | ref. 1000285     |
| · 1 time counter           | ref. 30880081    |
| · 1 active charcoal powder | ref. 121237.1609 |

## Scientific instrumentation

## Analytical techniques

## Chemical reagents

## Laboratory equipment

## Reducing matters and Previous inversion in:

Wines, Sparkling wines, Ciders  
Beers, Vinegars

...