



Sulphurous analysis kit GAB [1006000]

Ripper method for the detection of free and total sulphur in wines, musts, champagnes, vinegars, ciders, ...

Includes

- 1 x 500 mL starch [1006003]
- 1 x 1000 mL acid [1006005]
- 1 x 1000 mL alkali [1006002]
- 1 x 1000 mL valuation solution [1006004]
- Instructions manual

Characteristics

- Fast
- Easy method
- For red, rose, white wines, vinegars, ...
- Direct result

- Around 200 determinations
- Subjective result



- Scientific Instrumentation
- Analytical Techniques
- Laboratory Equipment
- Chemical Reagents

Free and total sulphur for red, rose, white wines, vinegars

...



Analytical Solutions