

Free and total SO₂ manual analyzer by Ripper method

Free and total SO₂ manual titrator by Ripper method able to detect titration's end point in red wines. Extremely easy and fast to use.

Characteristics

- Ripper method
- Compact design
- Simple and easy to use
- Very practical
- Without maintenance
- No need of calibration
- Power voltage: 12V
- Weight: 1,27Kg
- Dimensions: 16,5x17,5x13/50 cm



Equipment configuration

- 1 pyramid structure w/rod
- 1 burette clamp ref. 63204115
- 1 transformer ref. 0002065
- 1 beaker 100 mL p.p. ref. 409100
- 1 Toning tube ref. 3026056
- 1 stirrer magnet 6x20 mm ref. 137115
- 1 Toning tube cap ref. 47185050
- 1 starter reagents kit
- 1 instructions manual

Related products

- 2 graduated pipettes 5 mL ref. 3019013
- 1 graduated pipette 10 mL ref. 3019017
- 1 pipette 20 mL 2 m."A" ref. 3019023
- 1 burette 25 mL PTFE ref. 3019068

Scientific instrumentation

Analytical techniques

Chemical reagents

Laboratory equipment

Free and total SO₂ titrator

Wines, ciders, sparklings, ...