

## Free and total SO<sub>2</sub> manual analyzer by Ripper method

Manual titrator of free and total SO<sub>2</sub> by Ripper method able to detect titration's end point in red wines. Extremely easy to use.

### Characteristics

- Ripper method
- Compact design
- Easy to use
- Very practical
- Without maintenance
- No need of calibration
- Special modena vinegars
- Power voltage: 12V
- Weight: 2,15Kg
- Dimensions: 16,5x17,5x15/50 cm.



### Equipment configuration

- 1 pyramid structure with rod
- 4 glass tubes for toning ref. 3026056
- 4 rubber stoppers n°5 ref. 47185050
- 4 stirrer magnets ref. 137115
- 1 nut clip for burette ref. 63204115
- 1 inox rack for 12 tubes ref. R292
- 1 power voltage feeder ref. 0002065
- 1 beaker 100 ml p.p. ref. 409100
- 1 starter reagents
- 1 instructions manual

### Required material

- 2 graduated pipettes 5 mL ref. 3019013
- 1 graduated pipette 10 mL ref. 3019017
- 1 pipette 25 mL 2 m. "A" ref. 3019026
- 1 burette 25 mL PTFE ref. 3019068

Scientific instrumentation

Analytical techniques

Chemical reagents

Laboratory equipment

Wines

Sparkling wines

Ciders

Beers

Wines & Beers without alcohol

Vinegars

...